

**Labomiz**



# Class B Table-Top **Autoclave**



Email : [Info@Labomiz.com](mailto:Info@Labomiz.com) | Website : [www.Labomiz.com](http://www.Labomiz.com)

Class B-table top autoclave is a microcomputer controlled unit. It comprises of a three times better pre-vacuum drying function than the standard units. The built-in steam generator provides rapid cycles of saturated steam. Characterized with an automatic cooling air and steam exhaustion function at the end of the sterilizing process for user ease.

**Features:**

- ▶ Digital display with touch screen panel for ease of operation
- ▶ Built-in malfunction detector (over temperature, over voltage and low water level)
- ▶ Equipped with standard test interface
- ▶ Characterized with BD test, HELIX test and vacuum test
- ▶ Buzzing alarm indication and automated process end on completion of sterilization
- ▶ High quality SS liner with sealed covers
- ▶ Interlocking system for prevention of water spillage from tank

**Applications:**

Used in sterilizing dental tools, veterinary medicines, biological products, utensils, medical appliances, culture media, etc. in the fields of dentistry, biology, biochemistry, forensic, pharmaceuticals, cytology, immunology, genetics, bioengineering, health and sanitary, etc.

**Specifications :**

Model No.	MBTTA-1A	MBTTA-1B	MBTTA-1C	MBTTA-1D	MBTTA-1E
Capacity	8 L	12 L	16 L	18 L	23 L
Sterilization room size	Φ 170x360 mm	Φ 200x385 mm	Φ 230x390 mm	Φ 250x370 mm	Φ 250x470 mm
Exterior dimensions	370x345x565 mm	445x400x690 mm		480x465x690 mm	
Working Temperature	134 °C				
Maximum pressure	0.22 MPa				
Power	1500 w	1800 w		2200 w	
Voltage	110 V; 60 Hz				
	220 V; 50 Hz				
Weight	35 kg	46 kg	47 kg	52 kg	58 kg

# Labomiz Scientific

Email : [Info@Labomiz.com](mailto:Info@Labomiz.com) | Website : [www.Labomiz.com](http://www.Labomiz.com)